




Yamato Scientific
America



Parsortix[®] Platform


ISOLATE | ANALYZE | PERSONALIZE



cellbxhealth.com

Powerful cancer insights from blood

Parsortix® platform for Circulating Tumor Cell (CTC) enrichment



- Semi-automated, microfluidic cell enrichment platform
- Cell size and deformability-based cell capture
- Recover intact cells enabling multi-omic analysis using a wide range of molecular and imaging techniques
- Enables real-time biomarker assessment with the potential to support drug discovery and development, therapy selection and disease monitoring



[CLICK TO WATCH...](#)

Animation showing operation of the Parsortix cassette



The Parsortix platform uses patented microfluidic technology in the form of a single use cassette to enrich and recover viable CTCs and CTC clusters from blood, based on their less deformable nature and larger size.

This enables the isolation of a broad range of phenotypes including more invasive and treatment resistant mesenchymal cells and those undergoing epithelial-to-mesenchymal transition (EMT).



[CLICK TO WATCH...](#)

Click to watch patient blood flowing in the Parsortix cassette



The Parsortix workflow



1 BLOOD COLLECTION:
Collection of blood sample into vacutainer.



2 AUTOMATED BLOOD PROCESSING:
Prime the cassette. Insert your sample and press start – no preprocessing required. Blood is pumped through the cassette with minimal user input.



3 CELL CAPTURE:
Proprietary single use cassette captures CTCs within the 6.5 µm critical gap.



4 CELL IMAGING OR RECOVERY:
Examine cells in the cassette or reverse the flow to recover the cells in buffer for multiple downstream analysis techniques.

Cells recovered from the Parsortix platform are suitable for multi-omic analysis and are compatible with numerous downstream analysis techniques



Proteomic

Downstream applications such as FISH and IF imaging enable the assessment of specific nucleic acid targets and protein expression.

Transcriptomic

As intact cells, Parsortix-enriched CTCs are suitable for molecular applications including RNA sequencing.

Genomic



Parsortix-enriched CTCs are suitable for both targeted and whole-genome analyses, including PCR and NGS.

Functional

As viable cells, Parsortix-enriched CTCs can be cultured and used for preclinical studies e.g. xenografts.

Efficient performance...

- Capture individual CTCs and CTC clusters of varying sizes, independent of cell phenotype (*Figures 1 and 2*).
- High capture and recovery rates across multiple cancer types, with reliable cell enrichment from samples containing low cell numbers (*Figures 1 and 3*).



▲ *Figure 2*. Breast cancer CTC cluster containing epithelial, mesenchymal and epithelial-to-mesenchymal transition (EMT) phenotypes. Nuclear (blue), epithelial (green), mesenchymal (red), EMT (yellow/orange), blood cell (white).

▲ *Figure 1*. Efficient cell recovery of cell lines with different phenotypes. Strong linear correlation between recovery efficiency and low cell inputs between 2 and 100 cells per 7.5 mL blood.

CANCER TYPE	CELL LINE	CAPTURE EFFICIENCY %	RELEASE EFFICIENCY %
BREAST	SKBR3	94	80
	MCF7	74	73
	Hs578T	96	80
	BT549	64	91
LUNG	A549	72	90
PROSTATE	DU145	66	91
OVARIAN	CaOV3	84	84

▲ *Figure 3*. Efficient cell capture and release for multiple cell lines.


...with widespread adoption

Over 120 peer reviewed journal articles from independent study centers around the globe

- Most published marker independent CTC enrichment platform
- Parsortix platform utilized by 49 independent study centers and institutes globally
- Published literature investigates 19 cancer types representing >65% of all solid tumors
- Parsortix-enriched CTCs are suitable for a variety of multi-omic analysis techniques

VARIETY OF DOWNSTREAM ANALYSIS TECHNIQUES

RT-qPCR
 dPCR and ddPCR
 RNA-seq
 Immunofluorescence
 Fluorescence in situ hybridization (FISH)
 NGS/TGS
 WGA, WES and WTA
 Mass spectrometry
 Cell culture and xenotransplantation



CLINICALLY RELEVANT BIOMARKERS INCLUDING

EGFR	TP53
ESR1	AR
BRAF	AR-V7
KRAS	PIK3CA
PD-L1	DLL3
HER2	PSMA

NUMBER OF PUBLICATIONS BY CANCER TYPE: TOP 9

Breast	43
Lung	32
Prostate	18
Melanoma	8
Ovarian	7
Pancreatic	7
Colorectal	6
Brain	6
Head and Neck	6

Published use cases for Parsortix-enriched CTCs

- Prognostic assessment^{1,2}
- Biomarker assessment for personalized medicine³
- Minimally invasive, repeatable biomarker assessment to address tumor evolution⁴
- Combined assessment of ctDNA and CTCs for complementary biomarker assessment⁵
- Drug discovery and development⁶
- Unravelling the biology of cancer and the metastatic cascade⁷

Independent expert consensus highlights the Parsortix platform as a next-generation CTC technology for future clinical applications⁸

To discuss how we can support you, contact us at:

sales@cellbxhealth.com or visit cellbxhealth.com

CellBxHealth Europe Limited

10 Nugent Road,
Surrey Research Park,
Guildford, GU2 7AF,
United Kingdom
+44 (0)1483 343 434

customerservice@yamato-usa.com or visit
yamato-usa.com

Yamato Scientific America Inc.

2040 Corporate Court,
San Jose, California, 95131
United States
+1 408 235 7725

References

1. Davies, C. R. et al. *The potential of using circulating tumor cells and their gene expression to predict docetaxel response in metastatic prostate cancer.* *Front. Oncol.* 12, 1060864 (2023).
2. Obermayr, E. et al. *Gene expression markers in peripheral blood and outcome in patients with platinum resistant ovarian cancer: A study of the European GANNET53 consortium.* *Int. J. Cancer* *ijc.34978* (2024)
3. Ring, A. et al. *Circulating Tumor Cell Transcriptomics as Biopsy Surrogates in Metastatic Breast Cancer.* *Ann. Surg. Oncol.* 29, 2882–2894 (2022).
4. Bayou, N. et al. *Quantitative HER2 profiling on circulating tumor cells using an EpCAM-independent platform in metastatic breast cancer.* *Cancer Cell Int.* 25, 439 (2025).
5. Markou, A. N. et al. *Preoperative Mutational Analysis of Circulating Tumor Cells (CTCs) and Plasma-cfDNA Provides Complementary Information for Early Prediction of Relapse: A Pilot Study in Early-Stage Non-Small Cell Lung Cancer.* *Cancers* 15, 1877 (2023).
6. Kurzeder, C. et al. *Digoxin for reduction of circulating tumor cell cluster size in metastatic breast cancer: a proof-of-concept trial.* *Nat. Med.* 31, 1120–1124 (2025).
7. Bowley, T. Y. et al. *A Melanoma Brain Metastasis CTC Signature and CTC:B-cell Clusters Associate with Secondary Liver Metastasis: A Melanoma Brain-Liver Metastasis Axis.* *Cancer Research Communications*, 5(2), 295–308 (2025).
8. Nicolò, E. et al. *International expert consensus on the clinical integration of circulating tumor cells in solid tumors.* *European Journal of Cancer*, 231, 116050. (2025).



cellbxhealth.com



Yamato Scientific
America

yamato-usa.com