


# MARS<sup>®</sup> Bar Platform

VERSATILE PLATFORM SUPPORTING  
FLEXIBLE AND COST-EFFECTIVE  
CELL THERAPY DEVELOPMENT



**FAST AND  
FLEXIBLE**



**HIGH  
PERFORMANCE**



**COST-EFFECTIVE**





COST-EFFECTIVE

# Cost-Effective Cell Therapy Development from Lab bench to Clean room

## Cost-Effective

Column-free design reduces sample and reagent costs, very competitive running cost vs other market solutions

## Regulatory Ready

Designed for GMP environments with compliance-ready software and documentation

## GMP reagents

- CD4/CD8 nano biodegradable beads for T cell selection
- CD14 nano biodegradable beads for monocyte selection
- CD3/CD28 micro beads T cells selection and activation

## RUO

Fast and flexible protocol optimization


RUO reagents for cost-effective workflow

Up to 3 samples in parallel



# Advanced Magnetic Selection

Column-free cell selection  
with 50nm - 4.5 µm magnetic beads



## GMP Ready

21CFR Part 11 compliant with Audit Trail

Customizable and lockable protocols

Single use sterile fluids


ISO13485

FDA DMF submitted

CE, UL, RoHS, REACH



GMP Reagents





FAST AND FLEXIBLE

# Increased Throughput, Reduced Costs and Maximized Value



| Specifications               | Full LP       | ½ LP          | ¼ LP          |
|------------------------------|---------------|---------------|---------------|
| Pre-Label Volume             | 50 mL         | 25 mL         | 25 mL         |
| Label Incubation Time        | 30 min        |               |               |
| Pre-MARS® Volume             | 130 mL        | 65 mL         | 35mL          |
| <b>T Cell Selection Time</b> | <b>81 min</b> | <b>43 min</b> | <b>24 min</b> |
| Total MARS® Time             | 92 min        | 54 min        | 35 min        |



# END-TO-END PROCESS

## STEP

1

T cell selection reagents directly into LP sample



DIRECTLY FROM LEUKOPAK

OR



SPINNING

### Fast

Quick transition from unprocessed input sample to T cells in media without buffer exchange

### Customizable

Adaptable for quarter to full volume LP processing, maximizing the value of each run

### Reproducible

Consistent, high-quality T cell production with minimal handling

### Streamlined

From leukopak to CAR-T instreamlined workflow

## STEP

2

T cell selection  
T cells in complete media for G-Rex input



## STEP

3

Activation  
Expansion  
Fill & Finish

**OPTION 1**  
Cells eluted in medium, ready for immediate processing in a bioreactor



**OPTION 2**  
Cells eluted in buffer, ready for next step including freezing



### Hands-Free

Optimized for minimal hands-on operation

### Flexible

Compatible with traditional and fast CAR-T protocols for efficient, robust results

Experience a faster, more efficient,  
and cost-effective solution – MARS® delivers it all.



HIGH PERFORMANCE

# From CAR-T to Transformative Cell Therapies

## High Flexibility

Supporting a variety of sample types, volumes, and reagents

## High Purity, Recovery, and Viability


Top performance achieved for fresh or frozen leukopak samples, and directly from peripheral blood

## Direct Elution Options


Cell elution directly into culture media or buffers, achieving volume reduction and buffer exchange for immediate downstream processes

The new MARS® workflow enables the elution of immunomagnetically captured cells directly into the optimal medium, providing **rapid, centrifuge-free separation** and preparing cells for **immediate processing** (activation, expansion, etc) directly from MARS®. It also allows the elution of immunomagnetically captured cells in buffer, providing **rapid process** preparing cells for **freezing**.

### CD3+ CELL ISOLATION IN MEDIUM



### DEPLETION OF UNWANTED CELLS



Our proprietary, column-free technology delivers high purity and recovery of target cells, while effectively eliminating unwanted ones. This streamlined process using proprietary MARS® In-Situ Rinse saves time and **improves purity** and **recovery**, reducing manipulation steps and simplifying the process, **minimizing handling** and **saving time**.

# SPECIFICATIONS

## MARS® Bar Flex (Research use)

## MARS® Bar (Cell therapy development)

| SAMPLE                                   |                                                              |                                    |
|------------------------------------------|--------------------------------------------------------------|------------------------------------|
| Magnetic Cell labeling                   | ✓                                                            | ✓                                  |
| Containment                              | 5 mL, 15 mL, 50 mL tubes                                     | Bags                               |
| Sample processing                        | 1 sample each module Max 3 samples in parallel               | 1 sample by 3 parallel channels    |
| Sample types                             |                                                              |                                    |
| • Whole blood                            | ✓                                                            | ✓                                  |
| • Apheresis                              | ✓                                                            | ✓                                  |
| • Leukopaks                              | ✓                                                            | ✓                                  |
| • Frozen PBMC's                          | ✓                                                            | ✓                                  |
| • Bone marrow                            | ✓                                                            | ✓                                  |
| • Dissociated Tissue                     | ✓                                                            |                                    |
| REAGENTS & CONSUMABLES                   |                                                              |                                    |
| Isolation buffer                         | MARS® MAG buffer                                             |                                    |
| Isolation reagents                       |                                                              |                                    |
| • MARS® MAG lines (RUO)                  | ✓                                                            |                                    |
| • MARS® Ingenuity Line (RUO, GMP)        | ✓                                                            | ✓                                  |
| Fluidics                                 | Open-end tubing sets<br>Cleaning and sterilization protocols | Closed tubing set (Gamma radiated) |
| CELL ISOLATION                           |                                                              |                                    |
| Positive isolation                       | ✓                                                            | ✓                                  |
| • Direct from Whole Blood & Leukopak     | ✓                                                            | ✓                                  |
| Depletion                                | ✓                                                            | ✓                                  |
| Positive & negative tubing sets          | Same (program enabled change)                                |                                    |
| OPERATIONS                               |                                                              |                                    |
| Speed                                    | Protocol dependent; 0.5-6 mL/min                             |                                    |
| Column-free MARS® MAG in-flow technology | ✓                                                            | ✓                                  |
| Separation channels                      | Flex-BIBO scalable                                           |                                    |
| Redundancy                               | 3x Modules                                                   |                                    |
| Time to assemble tubing set              | <5 min*                                                      | <15 min*                           |
| Time to initiate isolation               | <2 min                                                       | < 8 min                            |
| Typical time to process 1e9 cells        | <15 min (75e6/mL)                                            |                                    |
| Capacity                                 | 0.5 - 45 mL per module                                       | 20 mL - 1L   Expandable >1L        |
| In-situ rinse for higher purity          | ✓                                                            | ✓                                  |
| Max total cells processable              | No practical limit                                           |                                    |
| Batched isolation                        | ✓                                                            | ✓                                  |
| Operation in bio-safety hood             | ✓                                                            | N/A                                |
| Additional configuration                 | Serial program**                                             | N/A                                |
| SOFTWARE                                 |                                                              |                                    |
| Pre-programmed protocols                 | ✓                                                            | ✓                                  |
| Adjustable & lockable parameters         | ✓                                                            | ✓                                  |
| Tiered user rights, Logged UI events     | ✓                                                            | ✓                                  |
| 21CFR Part 11 Compliant                  | ✓                                                            | ✓                                  |
| Encrypted logs                           | ✓                                                            | ✓                                  |
| INSTRUMENT                               |                                                              |                                    |
| Dimensions                               | 20.5" W x 16.5" D x 19.75" H                                 | 20.5" W x 16.5" D x 28" H          |
|                                          | 52cm W x 41 cm D x 50 cm H                                   | 52cm W x 42cm D x 74cm H           |
| Weight                                   | 62 lb / 28 kg                                                | 59 lb / 27 kg                      |

\* With standard training \*\* Customizable on demand

# Contact us

North America & International  
Applied Cells HQ

3350 Scott Blvd Bldg 6  
Santa Clara, CA 95054, USA  
Tel. 1-800-960-3004 EXT 1

[www.appliedcells.com](http://www.appliedcells.com) | [info@appliedcells.com](mailto:info@appliedcells.com)

## SINGAPORE

### SciMed Pte Ltd

Tel: (65) 6779 3388  
[sales@scimed.com.sg](mailto:sales@scimed.com.sg)  
[www.scimed.com.sg](http://www.scimed.com.sg)

## ITALY

### Euroclone S.p.A.

Tel. +39 02 38.19.51  
[info@euroclone.it](mailto:info@euroclone.it)  
[www.euroclone.it](http://www.euroclone.it)

## JAPAN

### Biomedica Solutions Inc

Tel 072-641-8140  
[info@bio-ms.com](mailto:info@bio-ms.com)  
[www.bio-ms.com](http://www.bio-ms.com)

## SPAIN AND PORTUGAL

### Izasa Scientific

Tel: +34 900 810 061  
[izasa@izasascientific.com](mailto:izasa@izasascientific.com)  
[www.izasascientific.com/es](http://www.izasascientific.com/es)

## INDIA

### Biotron

Tel: +91-22-6140 6400  
[info@biotronhealthcare.com](mailto:info@biotronhealthcare.com)  
[www.biotronhealthcare.com](http://www.biotronhealthcare.com)

## SOUTH KOREA

### DAON Biosciences, Inc

Tel +82 2 575 6227  
[www.daonbs.com/](http://www.daonbs.com/)

CZECH REPUBLIC, POLAND,  
SLOVAKIA, HUNGARY,  
ROMANIA

### Accela s.r.o.

Tel: +420 210 323 421  
[accela@accela.eu](mailto:accela@accela.eu)  
[www.accela.eu](http://www.accela.eu)

SWITZERLAND  
AND LIECHTENSTEIN

### Bucher Biotec AG

Tel. +41 (0)61 269 1111  
[info@bucher.ch](mailto:info@bucher.ch)  
[www.bucher.ch](http://www.bucher.ch)

## CHINA

### Nanjing Applied Cells Technology Co., Ltd.

Tel. +85 15032370112  
[Sales@appliedcells.com](mailto:Sales@appliedcells.com)  
[www.appliedcells.com](http://www.appliedcells.com)

## TURKEY

### CellTech

Tel: +90 322 232 47 27  
[info@celltech.com.tr](mailto:info@celltech.com.tr)  
[www.celltech.com.tr](http://www.celltech.com.tr)

## MEXICO

### RBC Química S.A. de C.V.

Tel: +52 55 7158 1833  
[ventas@rbcquimica.com](mailto:ventas@rbcquimica.com)  
[rbcquimica.com](http://rbcquimica.com)



ISO 13485:2016



## Learn more

<https://appliedcells.com/cell-therapy/>

For research use only. Not for use in therapeutic or diagnostic procedures. The MARS® Bar instrument and tubing set are designed, manufactured and tested under quality system certified to ISO 13485:2016. Not a medical device.

© Copyright 2025. All rights reserved. Applied Cells and MARS® are registered trademarks of Applied Cells, Inc. All other trademarks are the property of their respective owners.



B019E