

Stands


Stands and suitable oven models

Product code	Model	Suitable oven models
211180	ON30	DKL301C/311C, DKM300C/310C, DKN302C/312C, DNF301
211856	ON61	DX400/600 Series, DVS400/600 Series, DG400 Series, DKL400/600 Series, DKM400/600 Series, DKN400/600 Series, DNE400/600 Series, DNF400/600 Series
212348	OT42	DNE401/411, DNF401/411
212349	OT62	DNE601/611, DNF601/611
212477	OH41	DN411IE
212478	OH61	DN611IE
212801	ONS30	DX302C/312C
212802	ONS60	DX402C/412C/602C/612C, DR200/201
415464	OP43	DF411/412, DH411/412 (stand without caster)
415465	OP63	DF611/612, DH611/612 (stand without caster)
415466	OP46	DF411/412, DH411/412 (stand with caster wheels and stopper in front)
415467	OP66	DF611/612, DH611/612 (stand with caster wheels and stopper in front)

ON30



Dimension (mm)




ON30

OT42/62



Dimension (mm)




OT models

ON61



Dimension (mm)




ON61

OH41/61



Dimension (mm)




OH models

ONS60



Dimension (mm)




ONS60

OP43/63, OP46/66 (casters front side with two stoppers)



Dimension (mm) (Casters optional)



OP models

Stacking Kit


Stacking kit and suitable models

Product code	Model	Suitable oven models
212803	ODK80	DX302C/312C
212804	ODK82	DX402C/412C
212805	ODK84	DX602C/612C
212806	ODN26	400 Series of DNE, DNF (including models suitable for OD40)
212807	ODN28	600 Series of DNE, DNF (including models suitable for OD60)
212822	OD40	400 Series of DVS, DKL, DKM, DKN
212823	OD60	600 Series of DVS, DKL, DKM, DKN
213700	ODF48	DF412/612, DH412/612
281458	ODM44	DNF301

OD/ODN



Dimension (mm)




Overlaid bracket

ODK



Dimension (mm)



Overlaid bracket

Important Notes:

- The stacking units must belong to the same series (among 400 series or among 600 series)
- Do not stack 400 series on 600 series. If the upper 400 series is not fixed well, it is easy to topple
- Stacking of old models is forbidden. Lower unit must be new model.


Shelves

■ Stainless steel shelf & bracket and suitable oven models


Product code	Description	Model	Suitable oven models
211063	Shelf and bracket set stainless wire (loading up to 30 kg/shelf)	ODQ10	DF/DH412, DN411IE
211064	Shelf and bracket set stainless wire (loading up to 30 kg/shelf)	ODQ20	DF/DH612, DN611IE
211090	Shelf and bracket set stainless wire		DH650
211098	Shelf and bracket set stainless punch (loading up to 15kg/shelf)	ODQ30	DF/DH412, DN411IE
211099	Shelf and bracket set stainless punch (loading up to 15 kg/shelf)	ODQ40	DF/DH612, DN611IE
211854	Shelf and bracket set		DG800 Series
212068	Shelf and bracket set stainless punch		DKN/DX300/302, DNF301, DKM/DY300, DKL301/311
212095	Shelf and bracket set stainless punch		400 Series of DX, DVS, DKN
212192	Shelf		DP41/43/43C
212193	Shelf		DP61/63/63C
Q110204006	Shelf		DP83C
Q110204007	Shelf		DP103C/104C
212246	Shelf & bracket set stainless punch		400 Series of DVS, DG, DKL, DKM, DKN, DNE, DNF
212266	Shelf & bracket set stainless punch		600 Series of DX, DVS, DKL, DKM & 600 to 800 of DKN, DNE, DNF
212490	Shelf & bracket set stainless punch		900 Series of DKN, DNE, DNF
212678	Stainless steel wire shelf and bracket set. Withstand load up to 30 kg/shelf		DES830/DTS830
212679	Stainless steel punching shelf and bracket set. Withstand load up to 15 kg/shelf	ODE50	DES830/DTS830
212686	Stainless steel wire shelf and bracket set. Withstand load up to 30 kg/shelf		DE/DT411
212687	Stainless steel wire shelf and bracket set. Withstand load up to 30 kg/shelf		DE/DT611
212688	Stainless steel punching shelf and bracket set. Withstand load up to 15 kg/shelf		DE/DT411
212689	Stainless steel punching shelf and bracket set. Withstand load up to 15 kg/shelf		DE/DT611
212808	Shelf		DR200/201
212919	Basket-type shelf placed on top of std shelves. Withstand load up to 15kg/shelf	ODE12	DES830/DTS830
212924	Basket-type shelf placed on top of std shelves. Withstand load up to 15kg/shelf	ODT12	DF/DH412, DE/DT411
212925	Basket-type shelf placed on top of std shelves. Withstand load up to 15kg/shelf	ODT14	DF/DH612, DE/DT611
297071	Shelf		ADP21/200C/210C
297072	Shelf		ADP31/300C/310C
SHE-5680588	Tall shelf		SDP300/310
SHE-9751342	Short shelf		SDP300/310
SHE-5680563	Shelf		SDP400/410
SHE-5680562	Shelf		SDP610
YSA0000071	Shelf and bracket set		DF/DH832
YSA0000215	Shelf and bracket set		DF/DH1032




211063 / 211064




211090




211098 / 211099
212688 / 212689 / 212679




212068




212095




212192




212193



212246




212266



212490



212678 / 212686 / 212687



212808



212919 / 212924 / 212925



297071



297072

Earthquake Countermeasure

Seismic isolation rubber





Product code	296902
Material	Urethane elastomer
Max. load	~100 kg (by 4 pcs.)
Size	W50 x D50 x H5mm 1 set: 4pcs.

- Pasted at the bottom to prevent unit from falling
- Three layer structure made of urethane elastomer which absorbs 90% impact

Other Optional Accessories

Note: Listed optional accessories ARE NOT applicable to all oven models. Contact Customer Service for more details.

Cable port



Exhaust duct



Part no. R020603001
Butterfly flange
OD 100mm



Part no. R020603007
Double flange to suck in more air
OD 95mm x H 105mm



Exhaust duct for DF and DH
Part no. 213703 for DF / DH412
Part no. 213704 for DF / DH612

N₂ gas inlet device



External output terminal

