


Containers

Product code	Description	Dimensions	Suitable models
A990201305	Stainless rack	ø277 x H260mm	SK102C/112C
A990201304	Stainless rack	ø277 x H330mm	SK201C/211C
A990201303	Stainless rack	ø277 x H420mm	SK301C/311C
H060101047	Mesh basket	ø274 x H200mm	SN300C/310C/500C/510C
H060101048	Mesh basket	ø344 x H200mm	SQ500C/510C, SM520/530
H060103033	Mesh basket	ø344 x H300mm	SQ810C, SM820/830
H060601028	Mesh basket	ø424 x H200mm	SQL810C
H060602006	Mesh basket	ø424 x H300mm	SQL1010C
241085	Mesh basket (pitch 10 mm)	ø190 x H159mm	SM201/211
241086	Mesh basket (pitch 10 mm)	ø250 x H201mm	SM301/311/501/511
241087	Mesh basket (pitch 8.5 mm)	ø209 x H204mm	SM201/211, SN200C/210C
241088	Mesh basket (pitch 8.5 mm)	ø266 x H204mm	SM301/311/501/511
241089	Mesh basket (pitch 8.5 mm)	ø266 x H315mm	SM501/511
241090	Mesh basket with 2 stacking fittings	ø320 x H162mm	SQ500C/510C/810C, SM520/530/820/830
241091	Mesh basket with 3 stacking fittings	ø246 x H162mm	SM501/511, SN500C/510C
241092	Mesh basket with 2 stacking fittings	ø246 x H162mm	SM301/311, SN300C/310C
241093	Mesh basket with 2 stacking fittings	ø168 x H162mm	SN200C/210C, SM201/211
241094	Mesh basket with 2 perforated plates	ø270 x H590mm	SM501/511, SN500C/510C
241095	Mesh basket with 1 perforated plate	ø260 x H390mm	SM301/311, SN300C/310C
241096	Mesh basket with 1 perforated plate	ø200 x H390mm	SM201/211, SN200C/210C
241097	Mesh basket with 1 perforated plate	ø330 x H380mm	SQ500C/510C/810C, SM520/530/820/830
241083	Stainless solid basket	ø205 x H150mm	SM201/211, SN200C/210C
241084	Stainless solid basket	ø265 x H180mm	SM301/311/501/511, SN300C/310C/500C/510C
241150	Stainless bucket	ø208 x H203mm	SM201/211, SN200C/210C
241151	Stainless bucket	ø268x H203mm	SM301/311/501/511, SN300C/310C/500C/510C
241152	Stainless bucket	ø338 x H203mm	SQ500C/510C/810C, SM520/530/820/830

Mesh basket



H060103033

H060601028

241087

241088

241089

241090

241091

241092

241094

241095

241095

241096



H060101048

Stainless bucket



241150

241151

241152

Stainless rack



A990201303 / A990201304 / A990201305

Stainless solid basket




241083 / 241084


BOTTOM

Shelves


Product code	Punching shape	Suitable sterilizer models
212095	Round punch shelf & bracket set	SI401/402
212246	Round punch shelf & bracket set	SK401
212266	Round punch shelf & bracket set	SI601/602, SK601/801/811



212095



212246




212266

Stands


Product code	Stand models	Suitable sterilizer models
212802	ONS60	SI401/402/601/602
211856	ON61	SK401/601
212348	OT42	SK401
212349	OT62	SK601



ONS60



ON61



OT42/62

