





Aluminum block




Aluminum plate



Part Number	Code	Description	Weight	Compatible Containers
281432	OA167	Aluminum block	1.7 kg	Eggplant flask 200 mL
281436	OA171	Aluminum block	1.6 kg	Eggplant flask 300 mL
281566	OA172	Aluminum block	1.8 kg	Eggplant flask 500 mL
281567	OA173	Aluminum block	2.4 kg	Eggplant flask 1000 mL
281568	OA174	Aluminum block	2.4 kg	Eggplant flask 2000 mL
281572	OA175	Aluminum block	1.7 kg	Round flask 200 mL
281573	OA176	Aluminum block	1.7 kg	Round flask 300 mL
281574	OA177	Aluminum block	1.7 kg	Round flask 500 mL
281575	OA178	Aluminum block	2.4 kg	Round flask 1000 mL
281576	OA179	Aluminum block	2.4 kg	Round flask 2000 mL
281423	OA158	Aluminum plate	1.6 kg	Vial bottle Φ 12 mm 40 frame
281424	OA159	Aluminum plate	1.4 kg	Vial bottle Φ 15 mm 38 frame
281425	OA160	Aluminum plate	1.3 kg	Vial bottle Φ 17 mm 38 frame
281426	OA161	Aluminum plate	1.3 kg	Vial bottle Φ 18 mm 34 frame
281427	OA162	Aluminum plate	1.2 kg	Vial bottle Φ 21 mm 30 frame
281428	OA163	Aluminum plate	1.7 kg	Vial bottle Φ 30 mm 12 frame
281429	OA164	Aluminum plate	1.4 kg	Vial bottle Φ 35 mm 12 frame
281421	OA156	Base holder	0.4 kg	Adapter block 1 dress up
281422	OA157	Base holder	0.6 kg	Adapter block 3 dress up
281434	OA169	Adapter block	0.6 kg	Eggplant flask 10 mL
281433	OA168	Adapter block	0.6 kg	Eggplant flask 20 mL
281430	OA165	Adapter block	0.6 kg	Eggplant flask 30 mL
281431	OA166	Adapter block	0.6 kg	Eggplant flask 50 mL
281435	OA170	Adapter block	0.5 kg	Eggplant flask 100 mL
281577	OA180	Adapter plate	0.6 kg	Vial bottle Φ 12 mm 8 frame
281578	OA181	Adapter plate	0.6 kg	Vial bottle Φ 16 mm 6 frame
281579	OA182	Adapter plate	0.5 kg	Vial bottle Φ 24 mm 4 frame


Aluminum block
281568Aluminum plate
281424Base holder
281421 / 281422Adapter block
281430Adapter plate
281577

Part Number	Description	Code	Remarks
281395	Protective cover for MFD	OA154	Silicon protective cover protects the main unit from dirt and scattering of samples
281396	Protective cover for MFH	OA155	
281394	External temperature sensor	OA153	Sensor used when external temperature control is performed with MFH
281381	Pole set	OA143	Pole for fixing baths and mixers $\Phi 10 \times 480$ mm
281384	Container fall prevention frame	OA146	Variable frame to prevent containers such as beakers from falling when placed on the stirring table
281385	Stage for lab jack	OA147	Allows main unit to be placed on a lab jack
281587	Power cord (round terminal 2m)	OA183	Round terminal power cable
281382	Aluminum block handle	OA144	Handle to carry hot aluminum block with one hand
281383	Dual-handed aluminum block handle	OA145	Handles to carry hot aluminum blocks with both hands
231632	Muff	OLM44	$\Phi 5 \sim \Phi 13$ mm
231633	Muff	OLM46	$\Phi 6 \sim \Phi 17$ mm
231634	Muff	OLM48	$\Phi 9.5 \sim \Phi 29$ mm
231635	Double opening clamp	OLM50	Tightening adjustment range 3~55 mm, shaft $\Phi 10$ mm
231636	Double opening clamp	OLM52	Tightening adjustment range 3~80 mm, shaft $\Phi 20$ mm
222193	Glass tank for BOG100	OBO14	Mounted on MFH and used as an oil bath $\Phi 150$ mm 1.0L
222194	Glass tank for BOG200	OBO16	Mounted on MFH and used as an oil bath $\Phi 180$ mm 2.2L
281386	High magnetic agitator	OA148	Oval $\Phi 6 \times 15$
281390	High magnetic agitator	OA149	Octagon $\Phi 3 \times 13$
281391	High magnetic agitator	OA150	Octagon $\Phi 8 \times 13$
281392	High magnetic agitator	OA151	Octagon $\Phi 8 \times 38$
281393	Magnetic agitator	OA152	Micro $\Phi 2 \times 5$, set of 5
F-4028-02	Magnetic agitator	TB-20	$\Phi 7 \times 20$ 12 pcs.
F-4028-03	Magnetic agitator	TB-30	$\Phi 8 \times 30$ 12 pcs.
F-4028-04	Magnetic agitator	TB-40	$\Phi 8 \times 40$ 12 pcs.
F-4025-04	Magnetic agitator	A-43	$\Phi 13 \times 43$ 6 pcs.




Protective cover (MFH)

281396




Stage for lab jack

281385




Container fall prevention frame

281384




Dual-handed aluminum block handle

281383



Aluminum block handle

281382

High magnetic agitator
281390 / 281391 /
281392