

Desiccators

Designed to protect, store, or assist in transporting critical materials




860 SERIES

Dust and moisture free cabinets that safeguard valuable research specimens or sensitive electronic components. The status of components can be easily observed without opening the door.



860.5-SERIES

12" x 12" x 6"
304.8 x 304.8 x 152.4 mm



860-SERIES

12" x 12" x 12"
304.8 x 304.8 x 304.8 mm



861-SERIES

12" x 12" x 24"
304.8 x 304.8 x 609.6 mm



862-SERIES

18" x 18" x 18"
457.2 x 457.2 x 457.2 mm



863-SERIES

24" x 18" x 24"
609.6 x 457.2 x 609.6 mm

