









OVEN OVERVIEW

	Series	Model No.	Operating Temperature Range							Internal Capacity (L)	Program	Characteristics	
			0	100	200	300	400	500	600				700
Natural Convection	DX	302C/312C	RT+5~300°C							28	--	<ul style="list-style-type: none"> Economical High temperature Do not use fans. Heat rises by natural air convection for a slower heat flow 	
		402C/412C								74	--		
		602C/612C	RT+5~280°C							153	--		
	DVS	403C/413C	RT+10~260°C							99	Yes		
		603C/613C								162	Yes		
	DR	201	300~700°C							13.75	Yes		
	DG	400C/410C	RT+5~70°C							92	--		*DG840C/850C: Natural+Forced convection
		440C/450C								92	--		
		800C/810C								445	--		
		840C/850C*								445	--		
Forced Convection	DKL	301C/311C	RT+10~260°C							27	--	<ul style="list-style-type: none"> High level of air circulation, accuracy and uniformity Use fan motors for vertical air circulation providing a more uniform heat flow Built-in exhaust port Calibration offset function 	
		401C/411C								90	--		
		601C/611C								150	--		
	DKN	303	RT+10~260°C							27	Yes		
	DKV	400/410	RT+10~270°C							90	Yes		
		600/610								150	Yes		
		810								300	Yes		
		910	RT+10~220°C							545	Yes		
	DNE	401/411	RT+20~210°C							90	Yes		* DNF301/401/411/601/611 Two types of circulation: forced and natural convection
		601/611								150	Yes		
		811	RT+15~210°C							300	Yes		
		911								540	Yes		
	DNF	301	RT+15~260°C (Wind velocity: 1~10)							27	Yes		
		401/411								90	Yes		
		601/611								150	Yes		
811									300	Yes			
911									540	Yes			
Fine	DF	412	RT+15~260°C							91	Yes	<ul style="list-style-type: none"> Rapid & high volume of airflow Use forced convection for a horizontal air flow as opposed to vertical Very high uniformity, accuracy and performance Quick exhaust and cooling 	
		612								216	Yes		
		832	RT+15~200°C							512	Yes		
		1032								1000	Yes		
	DH	412	RT+15~360°C							91	Yes		
		612								216	Yes		
		650Z	RT+10~500°C							216	Yes		
		832	RT+15~300°C							512	Yes		
		1032								1000	Yes		
Vacuum	ADP	200C/210C	40~240°C							10	Yes	<ul style="list-style-type: none"> Handle sensitive samples at lower temperature Heat is evenly distributed from arrangement of the heaters against outer chamber walls Reduced oxidation 	
		300C/310C								27	Yes		
	SDP	300/310	RT+10~220°C							47.2	Yes		
		400/410								127.4	Yes		
		610	RT+15~220°C							264	Yes		
	DP	43C	40~200°C							91	Yes		
		63C								216	Yes		
		83C								512	Yes		
104C									1000	Yes			



OVEN OVERVIEW

	Series	Model No.	Operating Temperature Range							Internal Capacity (L)	Program	Characteristics
			0	100	200	300	400	500	600			
Inert	DN	411IE	RT+15~360°C							95	Yes	<ul style="list-style-type: none"> Creates non-oxidative environment Controllable nitrogen flow 
		611IE	RT+15~360°C							223	Yes	
Clean	DE	411	RT+30~260°C							91	Yes	<ul style="list-style-type: none"> Class 100 Adopts anti-fouling casters preventing wheel contamination during transportation Improved optional accessories and more customization options 
		611	RT+30~260°C							216	Yes	
	DT	411	RT+30~360°C							91	Yes	
		611	RT+30~360°C							216	Yes	
	DES	830	RT+30~260°C							327	Yes	
	DTS	830	RT+30~360°C							327	Yes	